

ORDER FORM

Prosthetic Socket

Patient ID: _____ **Side:** Left Right
Order Reference: _____ **Level:** TT TF¹ KD¹
Express: SY¹ HD/HP¹
Type: Check socket Definitive

Weight _____ kg **Acitivity level K1 – 4:** _____
Alignment^{2:} Standard³ CAD-file alignment⁴ CAD-file markings⁵
Extra material^{6:} Fibular head Tibial cut end
Distal adapter: Locking 65mm adult 4-hole⁷
 Locking 55mm small adult/pediatric 4-hole/TF valve⁸

Components: 30mm Pylon TI⁹ SS¹⁰
 30mm Tube Clamp TI⁹ SS¹⁰ 10°
 4-hole Male TI⁹ SS¹⁰ 10mm offset
 4-hole Slide/rotation Male
 Lamination adapter 72mm
 Locking 65mm Ball Bearing Water resist
 Locking 55mm small adult
 Locking 55mm small pediatric
Supplies: Loctite for plastic threads
 Scotchcast Plus 2" – 5cm¹¹

Notes:

1. Only available as check socket
2. Only applies to definitive socket with 4-hole adapter
3. Socket divided 50/50 proximally/distally in sagittal plane. 5° adduction in frontal plane. Adapter position determined accordingly
4. Adapter position according to alignment of CAD-file
5. Adapter position according to CAD-file markings
6. Extra material added on external socket wall for the

- possibility of adjustments
7. Prepared for suction/elevated vacuum
8. Prepared for suction/elevated vacuum. In addition TF-valve medially/distally or placed according to markings in CAD-file
9. Titanium
10. Stainless Steel
11. For securing of Check socket